



Description

Organic Acid Scale Remover is a concentrated, phosphate-free, organic acid detergent and descaler meant for cleaning and/or descaling cages, glassware, and mechanical equipment. **Organic Acid Scale Remover** is primarily designed for use inside washer disinfectors, cart washers, cage washers, glassware washers, and Clean-In-Place (CIP) applications. **Organic Acid Scale Remover** may also be used in chamber descaling applications.

Features

- High performance cleaning for Life Science
- Designed for use in immersion, pressure wash, and recirculating spray wash applications
- Effective against hard water scale, animal waste scale, beer stone, milk stone, etc.
- Low foaming and free rinsing
- Removes scale from cages, pans, cage washers, glassware washers, instruments, machines, and Clean-In-Place (CIP) systems
- Can also be used to passivate stainless steel

Directions for Use

Organic Acid Scale Remover should be dosed at 0.10-3.0 oz./gal. of water (0.8-24 mL/L). For descaling applications, **Organic Acid Scale Remover** should be dosed between 2.0-12.0 oz./gal. of water (16-96 mL/L) reverse osmosis (RO) or de-ionized (DI) water. Use temperatures can range from room temperature to 180°F (82°C). Always follow machine manufacturers' directions for use.

Safety

Be sure to don appropriate Personal Protection Equipment.

This product is intended for industrial and institutional use only. Not intended for consumer use.

Organic Acid Scale Remover is not compatible with some plastics and rubbers, and soft metals such as aluminum or brass.

Consult the Safety Data Sheet (SDS) for more detailed information pertaining to safety, handling, storage, and waste disposal methods.

DO NOT ADD BLEACH or any other chemical to **Organic Acid Scale Remover**. Combining chemicals outside of the directions for use could result in loss of efficacy or adverse chemical reactions.

The product is freeze/thaw stable. If frozen, thaw to room temperature and shake or stir.

Once the container is empty, rinse and dispose of properly. Do not reuse this container. Any detergent waste should be disposed of according to Federal, State, and/or local regulations.



Organic Acid Scale Remover

Organic Acid Glassware, Equipment,
and Cage Detergent and Descaler

Properties

Color: Colorless to light yellow liquid

Odor: Characteristic odor

Product pH: <4.0

Store below 120°F (49°C)

Use up to 35 months from date of manufacture

Package Sizes Available

5 Gallons (TH)

FICCB-CP-005T-8421

15 Gallons

FICCB-CP-0015-8421



Refer to Product Label and SDS for Additional Information.

314-771-6600 • 800-325-9218 • interconchemical.com

1100 Central Industrial Dr, St. Louis, MO 63110

2025 Intercon Chemical Company.

